










ISF1-Z / ISF1-ZC

ZC incubator shakers are optimized for cell cultivation

-  CO₂ control option available: essential for mammalian, plant cell cultures and algae
-  Heating and cooling
-  Controlled humidity option: essential for cultivation in microtiter plates or when cultivating in flasks for long periods
-  3 mm/1000 rpm option
-  Heated window and door frame with controlled humidity option
-  User-friendly operation: each parameter has its own control
-  Retrofitting possible
-  Touchscreen option available
-  **stack up to 3 shakers**
Easy to stack without the need for special tools or stacking kits



Automatic door and foot switch option available for easy handling, maximum comfort and automation



Stainless steel option





ISF-Z models and their features



Cooling / Heating
without CO₂ Control

Cooling / Heating
with CO₂ control



Basic





EcoDew®





Peltier



Peltier
with humidity control



Compressor



Compressor
with humidity control

