



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree Freezers

Sparkfree Freezer 124L RLVFO417

The Labcold Sparkfree range has a specially manufactured sparkfree interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it an important part of laboratory safety. For extra peace of mind this freezer is backed by a 2 year warranty.

For extra security a built in lock is fitted. The temperature display features high and low audio and visual temperature alarms along with a door open alarm.

WARNING: Do not store open containers of volatile substances in this freezer

Specifications

Product Code	RLVFO417
Gross Capacity	124 litres
Style	Undercounter
UK Warranty	2 years
External Dimensions (mm) HxWxD	830 x 600 x 600
Internal Dimensions (mm) HxWxD	680 x 460 x 410
Net Weight	45kg
Shelves	3

Specifications (cont.)

Temperature Range	-18°C to -25°C
Factory Set Point	-20°C
Door Style	Solid
Temperature Display	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Mains Failure Alarm	No
Self Closing Door	Yes
Access Porthole	Yes
Access Porthole Size	12mm
Manual Defrost	Yes
External Material	Powder Coated Steel
Internal Material	Powder Coated Steel
Lock Type	Key lock
Eco Friendly Refrigerant	Yes
Power Consumption Kwh/24hrs	0.8
Ambient Operating Temperature	16°C to 32°C

Features

- Sparkfree safety interior
- Door Lock
- Self closing door
- High/Low temperature alarm
- 3 shelves
- Manual defrost
- Access Port
- Temperature display
- Door Open alarm
- 2 Year parts and labour warranty (UK only)

Models

Model	Gross Capacity	Door Style	Lock Type
 Sparkfree Freezer 124L RLVFO417	124 litres	Solid	Key lock

Accessories



Labcold Stacking Kit