



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree Freezers

Sparkfree Freezer 416L RLCF1520

All Labcold chest freezers provide very stable storage conditions and slow temperature rise if power is interrupted.

They are rated for low energy consumption and include a lid lock for extra safety. Locating samples can be made simple and straight forward by use of a Labcold sample racking system – please ask us for a quote when you place your order.

Specifications

Product Code	RLCF1520
Gross Capacity	416 litres
Style	Chest
UK Warranty	1 year
External Dimensions (mm) HxWxD	865 x 1505 x 655
Internal Dimensions (mm) HxWxD	650 x 1300 x 430
Net Weight	67kg
Factory Set Point	-20°C




Specifications (cont.)

Temperature Range	-18°C to -25°C
Lid Style	Solid
Lock Type	Lid Lock
External Material	White Galv Steel
Internal Material	Steel
High Temperature Alarm	Visual
Baskets	2
Mains Indicator Light	Yes
Manual Defrost	Yes
Defrost Drain	Yes
Eco Friendly Refrigerant	Yes
Power Consumption Kwh/24hrs	0.72
Ambient Operating Temperature	16°C to 32°C

Features

- Sparkfree safety interior
- Mains indicator light
- Stable storage conditions
- Operator controlled defrost
- Defrost drain
- Lockable lid

Models

	Model	Gross Capacity	Lock Type
	Sparkfree Freezer 200L RLCF0720	200 litres	Lid lock
	Sparkfree Freezer 416L RLCF1520	416 litres	Lid Lock
	Sparkfree Freezer 556L RLCF2120	556 litres	Lid lock