



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree Freezers

Sparkfree Freezer 556L RLCF2120

All Labcold chest freezers provide very stable storage conditions and slow temperature rise if power is interrupted.

They are designed for low energy consumption and include a lid lock for extra safety.

This sparkfree chest freezer comes with baskets but Labcold can provide sample racking to help with organisation. Why not ask for details when you place your order?

WARNING: Do not store open containers of volatile substances in this freezer

Specifications

Product Code	RLCF2120
Gross Capacity	556 litres
Style	Chest
UK Warranty	1 year
External Dimensions (mm) HxWxD	865 x 1705 x 655
Internal Dimensions (mm) HxWxD	630 x 150 x 450
Net Weight	77kg
Factory Set Point	-20°C




Specifications (cont.)

Temperature Range	-18°C to -25°C
Lid Style	Solid
Lock Type	Lid lock
External Material	White Galv Steel
Internal Material	Steel
High Temperature Alarm	Visual
Baskets	2
Mains Indicator Light	Yes
Manual Defrost	Yes
Defrost Drain	Yes
Eco Friendly Refrigerant	Yes
Power Consumption Kwh/24hrs	0.80
Ambient Operating Temperature	16°C to 32°C

Features

- Stable storage conditions
- Defrost drain
- Lockable lid
- Operator controlled defrost
- Mains indicator light
- CFC and HCFC free

Models

	Model	Gross Capacity	Lock Type
	Sparkfree Freezer 200L RLCF0720	200 litres	Lid lock
	Sparkfree Freezer 416L RLCF1520	416 litres	Lid Lock
	Sparkfree Freezer 556L RLCF2120	556 litres	Lid lock