



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

### Sparkfree Freezers

## Sparkfree Freezer 66 litres RLVF0217

The Labcold Sparkfree range has a specially manufactured sparkfree interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it an important part of laboratory safety. For extra peace of mind this freezer is backed by a 2 year warranty.

### Specifications

Product Code	RLVF0217
Gross Capacity	66 litres
Style	Benchtop
UK Warranty	2 years
External Dimensions (mm) HxWxD	740 x 450 x 510
Internal Dimensions (mm) HxWxD	590 x 333 x 300
Net Weight	34kg
Shelves	3


## Specifications (cont.)

Temperature Range	-18°C to -25°C
Factory Set Point	-20°C
Door Style	Solid
Lock Type	Key lock
Temperature Display	Yes
High/Low Temperature Alarm	Yes
Self Closing Door	Yes
Door Open Alarm	Yes
Door Hinging	Right (can be reversed at time of manufacture)
Access Porthole	Yes
Access Porthole Size	12mm
External Material	Powder Coated Steel
Internal Material	Powder Coated Steel
Manual Defrost	Yes
Power Consumption Kwh/24hrs	1.4
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

## Features

- Sparkfree safety interior
- 3 shelves
- Temperature display
- Door open alarm
- 2 Year parts and labour warranty (UK only)
- Manual defrost
- Door Lock
- High/Low temperature alarm
- Self closing door

## Models

Model	Gross Capacity	Door Style	Lock Type
 Sparkfree Freezer 66 litres RLVF0217	66 litres	Solid	Key lock