



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

## Sparkfree Freezers

# Sparkfree Freezer RLVF1517

The Labcold Sparkfree range has a specially manufactured sparkfree interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it an important part of laboratory safety. For extra peace of mind this freezer is backed by a 2 year warranty.

For extra security a built in lock is fitted. The temperature display features high and low audio and visual temperature alarms along with a door open alarm.

## Specifications

|                                |                        |
|--------------------------------|------------------------|
| Product Code                   | RLVF1517               |
| Style                          | Freestanding           |
| Gross Capacity                 | 406 litres             |
| UK Warranty                    | 2 years                |
| External Dimensions (mm) HxWxD | 1875mm x 600mm x 700mm |
| Internal Dimensions (mm) HxWxD | 1690mm x 440mm x 500mm |
| Net Weight                     | 97kg                   |
| Shelves                        | 6                      |

## Specifications (cont.)

|                               |                |
|-------------------------------|----------------|
| Temperature Range             | -18°C to -25°C |
| Factory Set Point             | -20°C          |
| Door Style                    | Solid          |
| Lock Type                     | Key lock       |
| Temperature Display           | Yes            |
| Self Closing Door             | Yes            |
| High/Low Temperature Alarm    | Yes            |
| Door Open Alarm               | Yes            |
| Mains Failure Alarm           | No             |
| Access Porthole               | Yes            |
| Access Porthole Size          | 12mm           |
| Manual Defrost                | Yes            |
| Power Consumption Kwh/24hrs   | 3.16           |
| Eco Friendly Refrigerant      | Yes            |
| Ambient Operating Temperature | 16°C to 32°C   |

## Features

- Sparkfree safety interior
- Door Lock
- Self closing door
- High/Low temperature alarm
- 6 shelves
- Manual defrost
- Access Port
- Temperature display
- Door Open alarm
- 2 Year parts and labour warranty (UK only)

## Models

|                                                                                    | Model                      | Gross Capacity | Door Style | Lock Type |
|------------------------------------------------------------------------------------|----------------------------|----------------|------------|-----------|
|  | Sparkfree Freezer RLVF1517 | 406 litres     | Solid      | Key lock  |