



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree Fridges

Sparkfree Fridge 150L RLPRO517

The Labcold Sparkfree refrigerator range has a specially manufactured sparkfree interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it a significant aid to laboratory safety. For extra peace of mind it is backed by a 2 year parts and labour warranty.

For extra security an integral lock is fitted and the RLPRO517 features a NEW touch screen controller that display the temperature and features high and low audio and visual temperature alarms and a door open alarm.

WARNING: Do not store open containers of volatile substances in this fridge.

Specifications

Product Code	RLPRO517
Gross Capacity	150 litres
Style	Undercounter
UK Warranty	2 years
External Dimensions (mm) HxWxD	835 x 600 x 600
Internal Dimensions (mm) HxWxD	690 x 480 x 410
Temperature Range	0°C to 10°C
Factory Set Point	5°C

Specifications (cont.)

Internal Material	Powder Coated Steel
External Material	Powder Coated Steel
Net Weight	45Kg
Shelves	3
Lock Type	Key lock
Door Style	Solid
Temperature Display	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Access Porthole	Yes
Access Porthole Size	12mm
Auto Defrost	Yes
Power Consumption Kwh/24hrs	1.1
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

Features

- Sparkfree safety interior
- 3 adjustable shelves
- Temperature alarm and display
- Door Open alarm
- 2 year parts and labour warranty (UK only)
- Off Cycle Automatic defrosting
- Door Lock
- High/Low temperature alarm
- Access Port

Models

	Model	Gross Capacity	Lock Type	Door Style
	Sparkfree Fridge 150L RLPR0517	150 litres	Key lock	Solid

Accessories



Labcold Stacking Kit