



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree Fridges

Sparkfree Fridge 439L RLPR1517

The Labcold Sparkfree refrigerator range has a specially manufactured sparkfree interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it a significant aid to laboratory safety. For extra peace of mind it is backed by a 2 year parts and labour warranty.

For extra security an integral lock is fitted and the RLPR0517 features a NEW touch screen controller that display the temperature and features high and low audio and visual temperature alarms and a door open alarm.

WARNING: Do not store open containers of volatile substances in this fridge

Specifications

Product Code	RLPR1517
Gross Capacity	439 litres
Style	Free standing
UK Warranty	2 years
External Dimensions (mm) HxWxD	1875 x 600 x 700
Internal Dimensions (mm) HxWxD	1724 x 500 x 535
Net Weight	93Kg
Shelves	6

Specifications (cont.)

Temperature Range	0°C to 10°C
Factory Set Point	5°C
Door Style	Solid
Lock Type	Key lock
Temperature Display	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Rollers	4
Auto Defrost	Yes
Access Porthole	Yes
Access Porthole Size	12mm
Internal Material	Painted Aluminum
External Material	Painted Steel
Fan Assisted Cooling	Yes
Power Consumption Kwh/24hrs	2.55
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

Features

- Sparkfree safety interior
- 6 adjustable shelves
- Temperature alarm and display
- Door Open alarm
- Internal fan
- Off Cycle Automatic defrosting
- Door Lock
- High/Low temperature alarm
- Access Port
- 2 year parts and labour warranty (UK only)

Models

Model	Gross Capacity	Door Style	Lock Type
 Sparkfree Fridge 439L RLPR1517	439 litres	Solid	Key lock